

# INCUB-30 series

## Dry Bath Incubator



INCUB - 30S



INCUB - 30



INCUB - 30D - Dual block



Incubator	Power / Dimensions	Interface
Default Temperature: 37°C (±0.2°C) Default Time: 3 min / 5 min / 10 min / 15min Temperature Setting: 25°C to 60°C 3 Inbuilt Timers: T1 / T2 / T3 1 Location Open: Editable up to 99 hours / 59 min	Voltage: 110 - 220V / 50 - 60 Hz Power Supply: <60W Temperature: 10°C to 45°C Relative humidity: ≤80% Dimensions: 22cm X 27cm X 21cm	Screen: Alpha Numeric 20X2 lines Keyboard: Tactile Soft - Touch

Technical specifications subject to change without prior notice

Product manufactured by Micro Lab Instruments

An ISO 9001:2008 certified company